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Intended Purpose:

Extracorporeal blood circuit is intended to channel blood through external circuit for patients suffering from acute chronic renal failure and in hemodialysis.

A. Rinsing and air removal:

- (1) Unpack the Extracorporeal blood circuit.
- (2) Connect the arterial line to the physiological saline vial.
- *(3) Hang the drain bag which is at the end of venous line through correct place (hanging place of serum), and make sure the clamp is open.
- (4) Connect the arterial and venous line to the dialyzer and remove air completely by running physiological saline (recommended volume 1000 ml) at a flow rate of about 100 ml/min.
- * this one is applicable for blood line sets which have drain bag.

B. Priming:

Following rinsing and air removal, run about 500ml of physiological saline (eventually heparinized, 2000U/500ml) through the arterial line at flow rate about 100ml/min; then clamp the necessary sites on the arterial and venous lines. Complete all connections for dialysis. (Fig. 1)

C. Leakage test:

- (1) Prime the arterial and venous line completely with physiological saline by operating the blood pump; then stop the pump operation.
- (2) Clamp the arterial line near the dialyzer and the distal end of venous line with forceps.
- (3) Place the clamped distal end about 1 cm below the dialyzer and remove the forceps. (this results in application of a negative pressure of about 70mmHg to the blood compartment of the dialyzer.)
- (4) Examine for any or no continuous stream of bubbles in the venous header to check for leakage from the Extracorporeal blood circuit; if observed, replace the Extracorporeal blood circuit with a new one.

D. Start of dialysis:

- (1) Prepare the blood access site using sterile gloves. Insert the arterial and venous fistula needles with enough distance and in the suitable place. Note that the needles have sufficient distance from the fistula wall. Open the cap and clamp of the arterial line and after blood enters the tube, clamp the arterial needle. Do the same for venous line.
- (2) Connect the prepared blood access to the arterial line. Remove the forceps from the arterial and venous lines. While running dialysis at a flow about 500 ml/min, operate the blood pump at a flow rate of about 100 ml/min.
- (3) Make sure that no air bubbles remain in venous header or venous bloodline.
- *(4) Fully prime the arterial and venous lines including the dialyzer with blood by operating the blood pump, then stop the pump operation. Clamp the distal end of venous line and drain bag. Discard the bag in appropriate way.
- (5) Prepare the blood return site and connect to the venous line. After confirming that there are no bubbles in the line, remove the forceps from the line. After checking that there are no forceps on the lines and no line folding, operate the blood pump at a low flow rate. Take care not to apply excessive pressure to the lines and dialyzer to avoid leakage from the dialyzer and separation of each of the connection.
- (6) After confirming that there are no bubbles in the arterial and venous headers, turn the dialyzer 180° to allow removal of bubbles from dialysate. If bubbles are detected in the venous header before the turning, run blood at a prescribed flow rate for 5~10 minutes with the venous side kept upward.
- (7) Adjust the arterial drip chamber to about 1/3.
- (8) Adjust the venous drip chamber to about 2/3.
- *this sentence is practicable for blood line sets which have drain bags. For blood line sets which do not have a drain bag, clamp the distal end of venous line, then stop the pump operation.

E. Dialysis treatment-end and blood recovery

- (1) Stop the blood pump, clamp the arterial line and remove the line from the arterial blood access site; then connect the line to the physiological saline vial for blood recovery.
- (2) Unclamp the arterial line and run 100~200 ml of physiological saline at a flow rate of about 100ml/min to flush blood from the arterial and venous line and the dialyzer.

Doc.Code	Rev.No.:	Date Of Issue	Page
F72300RD05	R13	99/02/06 (25/Apr./2020)	1 of 4

(3) After blood recovery, discard the arterial and venous line and the dialyzer. Do not reuse them.

F. Indications

The Extracorporeal blood circuit is intended to facilitate channeling blood through external circuit outside the body when correctly installed to a Hollow Fiber Hemodialyzer, dialysate reservoir, and dialysis machine in patients suffering from acute or chronic renal failure for hemodialysis in patients.

G. Contraindications

Blood line set and fistula needle are sterilized with ETO, patients with sensitivity to ETO should not be exposed to Pharmed blood tubing.

H. Warning

- (1) Do not use for any others purpose than hemodialysis.
- (2) If foaming, blood leakage, blood coagulation and hemolysis occur during the use of this product, take appropriate measures according to a physician's instructions.
- (3) DO NOT reuse the device.
- (4) If the patient exhibits any abnormal symptoms such as discomfort, pruritus, urticaria, asthmatic reaction, hypertension and/or arrhythmia during the use of this product, take appropriate measures according to physician's instructions.
- (5) Pharmed medical industries will not be responsible for bad storage, bad transportation, abnormal usage and use of Extracorporeal blood circuit with any incompatible device with this circuit.
- (6) The expiry date refers to undamaged and properly stored product.
- (7) Kinked tubing can lead to hemolysis.
- (8) Laws restrict this device to sale by or on the order of a physician.
- (9) Pump segments can collapse partially at great pre-pump negative pressure or at reduced blood flow rates.
- (10) Discard the product, immediately after use by suitable means as contaminated medical waste and according to any prevailing environmental regulations.
- (11) Fistula Needle are DEHP Free and Bloodlines contains DEHP and it not recommended for children, pregnant and nursing women. therefore, PVC with DEHP plasticizer complies with the European Pharmacopoeia.
- (12) Sterilized with ETO and Non-pyrogenic

I. Caution

1. Caution before use

- a. Do not use if the package or if product is damaged.
- b. Use the product after opening the package immediately.
- c. Avoid air intake and contamination during rinsing and priming operations.
- d. check the seal of the connection and the seal of protection caps of unused inlets.
- e. check proper mounting of the pump tube and perform filling and rinsing of Extracorporeal blood circuit with a sterile pyrogen-free physiological solution.
- f. adhere these instruction for rinsing.
- g. make sure the cap of Fistula Needle is tightly closed.

2. Caution for storage

Store between (5°C~40°C) avoiding exposure to direct sunlight, severe vibration and high humidity.


3. It is recommended to heparinize the Extracorporeal blood circuit during the priming, the dosage of heparin and its mode of administration during dialysis are established by attending physician. The coagulation time should be checked regularly.

4. It is recommended to use of a dialysis machine equipped with a volumetric ultrafiltration unit, an accurate UF control system and a precise weigh for hemodialysis.

J. Connectors:

The connectors of Blood Line and Fistula Needle are female lure lock and they comply with the requirements of ISO 8637-1 and ISO 8637-2. The Blood Lines' connectors can connect to the blood ports of dialyzers.

Doc.Code	Rev.No.:	Date Of Issue	Page
F72300RD05	R13	99/02/06 (25/Apr./2020)	2 of 4

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
M. Manufacture

Name: Pharmed Medical Industries

Factory: 2nd Bustan, West Hafez Blvd, Eshtehard Industrial Zone, IRAN

Office: fourth floor, No. 44, Saadat Abad Ave., 1998898638, Tehran, IRAN

Email: Info@Pharmed.co

Website: www.Pharmaed.co

Telefax: +98 (26)377 77 620-5

Telefax: +98 (21) 220 73 102-3

N. European Representative

Name: NephroCan S.r.l.

Address: NephroCan S.r.l. VIA SAVOIA 78, 00198 ROMA, ITALY

Tel: +39 0685 2371 Fax: +39 06853 50187

www.NephroCan.com/EU-Rep