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
During cardiac bypass surgery

B. Indication

This set is called Hemoconcentrator, Hemoconcentration or Ultrafiltration and is used in bypass surgery if needed. The main target in this product is removing extra liquid and blood diluting.

C. Application

During cardiac bypass surgery, Intended for removing of extra water (used previous for blood dilution during priming oxygenator and other products in heart surgery or patient suffering unwanted extra liquid in blood circulation). This product is single use.

D. Side Effects

According to blood circulation outside of body, un- necessary utilizing, may result in blood hemolysis and related effects (releasing potassium in blood, potassium level in blood, exiting blood in the urine, broken blood cell sediment in kidney tubules and ...)

E. Product Segments

This product contains: Arterial Set, Venous Set, Hemoconcentrator and Drain Bag

Volume of blood pathway: Arterial line= 25 ml; Venous line= 25 ml

The Hemoconcentrator specification is as followed:

Hemoconcentration (Polyethersulfone 160HF) Idnt. Code: 4222350														
Type of Hemoconcentrator	UF Coefficient (Ultrafiltration) (ml/mmHg.h)	KoA (ml/min)	Clearances (ml/min)										Blood priming Volume(ml)	Surface Area (m ²)
			Urea		Creatinine		Phosphate		Vitamin B ₁₂		Inulin			
Polyethersulfone 160 HF	62	1167	$Q_B = 200$	$Q_B = 300$	$Q_B = 200$	$Q_B = 300$	$Q_B = 200$	$Q_B = 300$	$Q_B = 200$	$Q_B = 300$	$Q_B = 200$	$Q_B = 300$	89	1.6
Hemoconcentration (Polysulfone 160 HF) Idnt. Code: 4222350														
Polysulfone 160 HF	55	1145	195	270	191	252	183	233	142	165	108	120	86	1.6

Specification and performance data are estimated at Q_B : 200/300 ml/min, Q_B : 500 ml/min, Q_F : 0 ml/min, T: 37°C.

Performance data were measured in vitro according to standards BS EN ISO 8637-1 and BS EN ISO 8637-2

Typical values obtained with an individual batch of fiber, the clinical use may determine difference in results in relation to different ultrafiltration and measuring technique possible verification between batches of fibers.

UF measurement using Bovine/Human blood (Hct 32%, Protein 60 g/l)

The tolerance for clearance numbers is $\pm 10\%$

HF: High Flux

F. Contra-indications

This product must not be used if not needed or in patients having Hypovolemia.

G. Before Using

- Don't use if the product or packing is damaged.
- close the discharge valve of drain bag before use
- Avoid air intake and contamination during rinsing and priming operation
- Use immediately after priming.
- Knowing hemoglobin, hematocrite, height, weight and patient's temperature is necessary.
- knowledge about the serums and drugs used for blood circulation outside of body.

Doc. Code	Rev. No.:	Date Of Issue	Page
F72300RD52	R08	00/01/23(12/Apr./2021)	1 of 2

H.Lines and Cannulation

Arterial and venous cannulation are performed by surgeon and are fixed with sutures in their places. Cannulation is used for blood supplement outside of body with placing two cannulas in upper and lower vertical vein or placing a single cannula in tow steps inside of corridor. Arterial Cannulation is done by placing a special cannulas inside of ascending aorta and then connected to outside of body lines.

I.During Usage

- Control oxygenator, arterial and venous pressure continuously
- Check drain bag
- Ensure patient temperature

J.After Usage

The hemoconcentrator (dialyzer) and set should be disposed in a safe manner as means as contaminated medical waste and according to any prevailing environmental regulations.


K.Storage

This Product should be stored between +5 to +35°C. Avoid sunlight, intense radiation and high humidity.

L.Connections

The blood ports of hemofilter comply with the requirements of BS EN ISO 8637-1 and BS EN ISO 8637-2 .They can connect to female luer lock connectors of blood line set.

Arterial and venous sets are connected to Hemoconcentrator filter through connectors (Red and Blue) and connected to oxygenator in the other side which supplies oxygen for bloods outside of body and acts as an artificial lung. It should be noted that the blue connector should be connected in the side with drain bag. Oxygenator has a pump and set out in pressure $Q_B=200-300$ ml/min and 40-60 mmHg. This object is connected to oxygenator and moves with heart pump.

M.Manufacture

Name: Pharmed Medical Industries
 Factory: 2nd Bustan, West Hafez Blvd, Eshtehard Industrial Zone, IRAN Telefax: +98(26) 3777 7620-5
 Office: fourth floor, No. 44, Saadat Abad Ave., 1998898638, Tehran, IRAN Telefax: +98(21) 2207 3102-3
 Email: Info@Pharmed.co Website:www.pharmed.co

N.Symbols

	Manufacturer		Do not use if package is damaged		Max TMP		The side up		Read the instructions before use
	Fragile, handle with care		Keep away from sunlight		Caution		Do not open with sharp object		Single use
	Max stuck		Non - pyrogenic		Sterilized with ETO		Keep dry		Keep at 5°C~35°C

Note: Pharmed Medical Industries is not responsible for any shipping accidents, bad storage and abnormal use.

Doc. Code	Rev. No.:	Date Of Issue	Page
F72300RD52	R08	00/01/23(12/Apr./2021)	2 of 2