ICD GROUP

Instruction for use of RenaCart, RenaBag & RenaBag S (Sodium Bicarbonate Powder Cartridge and Bag)

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A. Intended Use

RenaCart (Sodium Bicarbonate Powder Cartridge), RenaBag and RenaBag S (Sodium Bicarbonate Powder Bag) are used to prepare hemodialysis solution with bicarbonate base, which with water and Concentrated Hemodialysis Acidic Solution making dialysis fluid. (According to the instructions of the dialysis machine)

RenaCart, RenaBag and RenaBag S are intended to be used during Hemodialysis and online Hemodiafiltration.

B. Cautions

1. The information in this manual must be read carefully and understood before using the product.
2. Do not refill the cartridge and bag of sodium bicarbonate and do not reuse them.
3. The fluid obtained from Sodium Bicarbonate Powder should be consumed immediately after preparation.
4. Do not use if the cap is missed or opened.
5. ISO-quality water meeting the requirements of BS EN ISO 23500-3 shall be used to dilute the concentrate to make dialysis fluid.

C. Warning

1. This product is single use.
2. This product is for one dialysis treatment session. Discard any unused residual powder or solution.
3. This product is non-pyrogenic and non-sterile.
4. The product variety table is in part M.
5. Storage conditions: Store at 5 to 30°C, Keep out of freezing.
6. Do not use it, if the bag or cartridge is damaged or suspected.
7. This product can only be used with Concentrated Hemodialysis Acidic Solution recommended in the dialysis machine manual.
8. Do not initiate dialysis without having followed the instructions the operator's manual of the corresponding dialysis machine regarding any recommended tests to verify the proper composition of the dialysis fluid.
9. Do not make unauthorized changes to the bag and cartridge of sodium bicarbonate or to the dialysis machine, which will lead to serious consequences on the therapeutic effects of hemodialysis or lead to improper functioning of the device.
10. The products shall be used only by trained nurses (personnel) and/or physician who are totally familiar with procedure and exactly according to the instruction for use.
11. Pharmed Medical Industries will not be responsible for bad storage, bad transportation, abnormal use and use of the product with any incompatible device.
12. The bag and cartridge must be disposed, after use, in a safe manner as mean as contaminated medical waste and according to the environmental health regulations.

D. Description

1. RenaCart (Sodium Bicarbonate Powder Cartridge) is made Polypropylene (PP) and contain 550g or 650g or 750g dry Sodium Bicarbonate Powder. RenaBag and RenaBag S (Sodium Bicarbonate Powder Bag) is made Polyamide/Polyethylene (PA/PE) and contain 550g or 650g or 750g dry Sodium Bicarbonate Powder.

2. When the cartridge or bag is placed on the special holder, purified water conduct by the dialysis machine through the product, producing a saturated solution of sodium bicarbonate. The dialysis machine mixes the sodium bicarbonate powder solution with purified water and the acid concentrate to produce a bicarbonate-based dialysis fluid.

3. Both the powder and liquid solution, obtained from sodium bicarbonate cartridges or bags, meet the requirements of the European and US pharmacopeias.

E. Contraindications

There are no absolute contraindications for bicarbonate dialysis.

Sodium Bicarbonate Cartridge or Bag must be used only on the directions of a physician who has evaluated all the pertinent features of the product in relation to individual patient.

F. Side effects

Some undesirable effects related to dialysis can occur, for instance hypotension, nausea, vomiting and cramps.

G. Precautions for use

1. RenaCart, RenaBag and RenaBag S (Sodium Bicarbonate Powder Cartridge and Bag) which are suitable for all acidic formulas as per the ration stated on the label of container of concentrated hemodialysis solution.

Attention: RenaCart, RenaBag & RenaBag S contains 550g or 650g or 750g of Sodium Bicarbonate Powder which is suitable to dissolve in 8 liters of purified water and sufficient for one dialysis treatment of 4 to 5h. At a dialysate flow rate of 500ml/min and sodium bicarbonate concentration of 35 mEq/l.

2. Acid and basic bicarbonate hemodialysis concentrates should only be diluted for short time before use.

3. The composition of the final dialysis fluid for use must be closely monitored.

4. Use just only for a dialysis session. The bicarbonate solution is highly susceptible to microbial contamination so discard any unused solution.

H. General instructions for use

See operator's manual of instructions for the corresponding hemodialysis machine.

I. Specific Instructions for installation

RenaCart (Sodium Bicarbonate Powder Cartridge)

1. Open the two latches of the cartridge holder.
2. Insert the lower cartridge connection into the hole of the lower latch.
3. Press down the upper latch and insert the upper cartridge connection into the hole of the upper latch.
4. Press down the upper latch until you hear a noise indicating that the cartridge is locked in place. Rotate the cartridge 90° to ensure the free flow of the liquid.

RenaBag and RenaBag S (Sodium Bicarbonate Powder Bag)

1. Remove the cap.
2. Hold vertically with port applied toward the machine.
3. Insert the connector of the bicarbonate bag to its appropriate nozzle of the dialysis machine.
4. Press forward till click that means complete fixation of the bag to the machine.
5. The machine will automatically detect the bicarbonate bag and will start dissolving the Sodium Bicarbonate Powder to form the bicarbonate solution needed.

Doc. Code	Rev. No.:	Date Of Issue	Page
F72300RD04	R13	01/07/18(10/Oct./2022)	1 of 2

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J. General instructions to disconnect: After disconnecting the patient, proceed as follows:

RenaCart (Sodium Bicarbonate Powder Cartridge)

1. If the control program of the machine includes the cartridge discharge procedure before the cartridge is removed from the holder after dialysis, perform this procedure first and then remove the cartridge from its holder. Otherwise, proceed as per the following points.
2. Lift the upper latch.
3. Place a finger on the upper cartridge opening, at the same time lift the cartridge out and twist it if appropriate. You can now remove the cartridge.
4. Clean the cartridge holder and close both latches.
5. Proceed further in accordance with the instructions of the directions for use of the corresponding hemodialysis machine.

RenaBag and RenaBag S (Sodium Bicarbonate Powder Bag)

1. After finishing the hemodialysis session, press on the bottom beside the nozzle of bicarbonate bag in the machine for RenaBag and for RenaBag S remove the cover of seat location of RenaBag S.
2. Emptying the bag from the remaining solution by pressing on the control panel of the software of hemodialysis machine (Empty bag) to remove all the residual solution in the bag before disconnecting the bag from the machine.
3. Remove the Bag.

K. Type of hemodialysis machine

RenaBag shall be used with Fresenius hemodialysis machine 4008 B and RenaBag S shall be used with Fresenius hemodialysis machines 4008 S & 5008 which are designed to use bicarbonate dialysis.

RenaCart is used in the following machines designed to use bicarbonate dialysis:

Dialysis Machine Type	Nipro	Bellco			Hospal		Gambro			B-Braun	
Model	Nipro Sordial 55 Plus	Bellco formula 2000	Bellco formula Therapy	Bellco formula Plus	Hospal Integra	Hospal Montreal	AK 95	AK 200	AK 96	B-Braun Dialog	B-Braun Dialog+

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O. Manufacture

Name: Pharmed Medical Industries
 Factory: 2nd Bustan, West Hafez Blvd., Eshtehard Industrial Zone, IRAN
 Telefax: +98 (26)3777620-5
 Office: fourth floor, No. 44, Saadat Abad Ave., 1998898638, Tehran, IRAN
 Telefax: +98 (21) 22073102-3
 Email: info@pharmed.co Website: www.pharmed.co
 Customer Care Services: +989123346849 +98 (21)22143509

P. European Representative

NephroCan
 NephroCan S.r.l., via savovia 78, 00198 ROMA, ITALY
 Tel Office: +39 0685 2371 Fax: +39 06853 50187
 www.NephroCan.com/EU-Rep

Doc. Code	Rev. No.:	Date Of Issue	Page
F72300RD04	R13	01/07/18(10/Oct./2022)	2 of 2