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**A. Intended Use**

**1. Acidic formulae; A-Concentrate**

Concentrated Hemodialysis solution Acidic Formula A-concentrate [dialyzing fluid will be mixed with physiological or higher concentrations of bicarbonate B-concentrate as buffer] For extracorporeal bicarbonate dialysis use A-concentrate only in combination with B-concentrate (NaHCO<sub>3</sub> 4.8% or 6.6%) in the given dilution or with a Bicarbonate container (RenaCart, RenaBag or Renabag S); is intended for removing undesirable substances (toxic substance such as Urea, Creatinine, phosphate, and other uremic toxins) out of the blood by exchange of solutes through a diffusion process, and separated from blood by a semi-permeable dialysis membrane (hollow fiber membrane) and including online procedures only in combination with HD monitor capable and intended for online filtration sufficient for dialysis treatment duration of 4-6 hours (see relevant dialysis machine operators manual).

**2. Acetate formulae; without buffering solution**

Concentrated Hemodialysis solution Acetate Formula [dialyzing fluid without bicarbonate, using acetate as a buffer] For extracorporeal bicarbonate dialysis use Acetate formulae only without buffering solution in the given dilution is intended for removing undesirable substances (toxic substance such as Urea, Creatinine, phosphate, and other uremic toxins) out of the blood by exchange of solutes through a diffusion process, and separated from blood by a semi-permeable dialysis membrane (hollow fiber membrane) and including online procedures only in combination with HD monitor capable and intended for online filtration sufficient for dialysis treatment duration of 4-6 hours (see relevant dialysis machine operators manual).

**B. Indication**

**1. Acidic formulae; A-Concentrate**

When mixed with B-Concentrate, it is intended to be used for dialysis patients who are suffering from Acute and chronic renal failure, hyperhydration, and intoxication, correction of acid/base metabolism, and electrolyte status.

Concentrated Hemodialysis Acidic Solution	Potassium Electrolyte Concentration (mmol/l)
Concentrated Hemodialysis Acidic Solution I Idnt.Code: 4316205	1
Concentrated Hemodialysis Acidic Solution II Idnt.Code: 4316210	2
Concentrated Hemodialysis Acidic Solution III Idnt.Code: 4316215	3

**2. Acetate formulae; without buffering solution**

When mixed with B-Concentrate, it is intended to be used for dialysis patients who are suffering from Acute and chronic renal failure, hyperhydration, and intoxication, correction of acid/base metabolism, and electrolyte status.

**C. Side Effects**

Use of Concentrated Hemodialysis Solutions (both formulae; Acidic and Acetate) may lead to some contraindications these are as follows: Hypotonia, nausea, vomiting and spasms have been reported.

**D. Contra-indications**

- That the solution is not to be used for intravenous infusion.
- Use of concentrated Hemodialysis solutions (both formulae; Acidic and Acetate) may lead to some contraindications these are as followed:
  - Acidic Hemodialysis Concentrate I don't use Hypokalemia (applicable for potassium-free concentrate), and unmanageable coagulation disorders.
  - Acidic Hemodialysis Concentrate III don't use Hyperkalemia (applicable for potassium containing concentrates).
  - For patients with haemodynamic instability, an alternative method of extracorporeal treatment may be indicated.

**E. Caution and Warning**

- Both formulae of concentrated Hemodialysis Solutions (A-concentrate with B-concentrate, and Acetate) must be diluted with Purified water immediately before use in accordance with the relevant national requirement. Discard residues.
- Always check the concentration of the ready to use dialysis solution. Use only clear solution.
- Do not use from damaged containers.
- Concentrates are Non-Sterile and non-pyrogenic.
- Keep at 4-30°C. (Short-term exposure to 40°C will not harm to acid concentrated).
- Keep out of light and freezing.
- Connect Red connection tube with A-Concentrate, and Blue connection tube with B-concentrate.
- In case of treatment interruptions (Multiple Hemodialysis sessions); make sure that the container of the concentrated Hemodialysis solution is covered well, and keep it clean and away from any source of contamination.
- For product performance; refer to the concentration of Hemodialysis solution mentioned on the product label.
- To be used only by prescription.

**F. Manufacture**

Name: Pharmed Medical Industries


Factory: 2<sup>nd</sup> Bustan, West Hafez Blvd, Eshtehard Industrial Zone, IRAN    Telefax: +98 (26) 3777 7620-5

Office: fourth floor, No. 44, Saadat Abad Ave., 1998898638, Tehran, IRAN    Telefax: +98 (21) 2207 3102-3

Email: info@pharmed.co    Website: www.pharmed.co

Customer Care Services: +989123346849    +98(21)22143509

**G. Symbols**

	Non sterile		Manufacturer		Max stuck		High-Density Polyethylene		Read the instructions before use
	Do not step		Fragile, handle with care		The side up		Keep dry		

Doc. Code	Rev. No.:	Date Of Issue	Page
F72300RD09	R04	1401/02/24(14/May/2022)	1 of 1