




Blood Bags

Pharmed Medical Industries
(Member of ICD Group)

F72300RD106 R01

PHARMed
Medical Industries

EN ISO 13485



PHARMED

Medical Industries

Blood Bags

Product Code	Configuration	Collection vol.	CPD vol.	Transfer Bag 1	Transfer Bag 2	Transfer Bag 3
4512240	Single Blood Bag	450 ml	63 ml	—	—	—
4512265	Double Blood Bag	450 ml	63 ml	Empty Transfer bag suitable for storage of plasma	—	—
4512230	Triple Blood Bag	450 ml	63 ml	Empty Transfer bag suitable for storage of plasma and 5 days platelet storage		—
4512272	Quadruple Blood Bag	450 ml	63 ml	Empty Transfer bag suitable for storage of red cells, plasma and 5 days platelet storage		

Product Code	Configuration	Collection vol.	CPDA-1 vol.	Transfer Bag 1	Transfer Bag 2	Transfer Bag 3
4512250	Single Blood Bag	450 ml	63 ml	—	—	—
4512260	Double Blood Bag	450 ml	63 ml	Empty Transfer bag suitable for storage of plasma	—	—
4512220	Triple Blood Bag	450 ml	63 ml	Empty Transfer bag suitable for storage of plasma and 5 days platelet storage		—
4512280	Quadruple Blood Bag	450 ml	63 ml	Empty Transfer bag suitable for storage of red cells, plasma and 5 days platelet storage		

Product Code	Configuration	Collection vol.	CPD vol.	SAGM vol.	Transfer Bag 1	Transfer Bag 2	Transfer Bag 3
4512210	Triple Blood Bag	450 ml	63 ml	100 ml	Empty Transfer bag suitable for storage plasma	SAGM Transfer bag suitable for storage red cells	—
4512290	Quadruple Blood Bag	450 ml	63 ml	100 ml	Empty Transfer bag suitable for storage of plasma	Empty Transfer bag suitable for storage 5 days platelet storage	SAGM Transfer bag suitable for storage red cells

Pharmed Medical Industries

4th Floor, No.44, Sa'adat Abad St., Tehran, Iran

Tel: +98 (21) 2207 3102

Fax: +98 (21) 2207 3103

